
LSMMed

DCD Cerclage

Product description
& surgical technique

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Indications

- Femoral fractures.
- Humeral and tibial fractures.
- Diaphyseal reconstruction in revision hip arthroplasty.
- Diaphyseal reconstruction in revision shoulder arthroplasty.

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Materials

Cable: titanium alloy Ti6Al4V in conformity with ASTM F136 standard.

Spacers and crimp: titanium type 2 in conformity with ISO 5832-2 standard.

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Technical features

- Thanks to the spacers, it is possible to avoid the damage to the periosteal blood flow.
- Easier bone healing.
- MRI compatible.
- Biocompatibility.
- Reduced allergic reactions (no Nickel into the alloy).
- Compatible with other titanium implants (plates, screws, revision stems, ...).

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Clinical cases

Casi clinici

Revision hip arthroplasty:
periprosthetic fracture treated by using
DCD cerclage.

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

Surgical technique

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1.

1. Position the patient according to the surgical approach and reduce fracture. Select the appropriate cable passer, HME 869-1 or 869-2, so that the instrument passes around the bone without causing significant damages to soft tissues or excessive stripping of the periosteum.

Place the cable passer around the bone. Insert a suture trial wire in the cable passer to determine the needed cerclage cable length and therefore spacers placement.

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
2.


2. Assemble the cerclage cable for the implant. Insert first the cable crimp all the way down up to the ball end of the cable and with the four points on the under surface of the crimp in contact with the instrument table. Then insert the spacers by distancing them according to the previously determined cable length and with the notch in parallel to the cable crimper pliers HME 867, as shown in the picture; tighten them on the cable by using the above-mentioned crimper. Number of spacers to use is variable (one package contains 6 spacers).

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Surgical technique

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3.

3. Insert the cable end with no spacers in the end hole of the cable passer, as shown in the picture.

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
4.

4. Remove the cable passer leaving the cable wrapped around the bone. Insert the end of the cable through the free hole of the crimp and tighten up to spacers and the four points on the under surface of the crimp are in contact with bone.

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
5.

5. Prepare the tension of the cerclage cable. Turn counterclockwise the fluted rear knob on the cable tensioner HME 866 up to the tension reference line disappears distally, as shown in the picture. Unscrew the lateral locking screw and insert the cable. Advance the tensioner bit up to the crimp and screw the lateral screw to lock the cable in the tensioner.

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Surgical technique

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
6.

6. Turn clockwise the fluted rear knob on the cable tensioner until the suitable tension is reached (tension reference is shown by the markings on the tensioner 10-40 kg); tension of cerclage cable should not exceed 40 kg to avoid possible bone injuries.

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7.

7. When cable tension is reached, secure the cerclage cable by placing the cable crimper HME 867 on the crimp and tightening it firmly. Unscrew then the lateral screw and remove the tensioner; cut the loose end of the cable by using the nipper HMG 0917.

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Bibliography

Berend KR, Lombardi AV Jr, Mallory TH, Chonko DJ, Dodds KL, Adams JB.
"Cerclage wires or cables for the management of intraoperative fracture associated with a cementless, tapered femoral prosthesis: results at 2 to 16 years."
J Arthroplasty. 2004 Oct;19(7 Suppl 2):17-21.

Sandhu R, Avramidis K, Johnson-Nurse C.
"Dall-Miles cable and plate fixation system in the treatment of periprosthetic femoral fractures: a review of 20 cases."
J Orthop Surg (Hong Kong). 2005 Dec;13(3):259-66

Jarit GJ, Sathappan SS, Panchal A, Strauss E, Di Cesare PE
"Fixation systems of greater trochanteric osteotomies: biomechanical and clinical outcomes"
J Am Acad Orthop Surg. 2007 Oct;15(10):614-24.

Info

LSM-Med Srl
reserves the right
to make changes.

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OVERMED Srl