RSB Implants

Resorbable implants
Impianti riassorbibili

Product description
& surgical technique

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Technical features

The application of resorbable materials for orthopaedic, general, maxillofacial surgery, etc. dates back to the beginning of 80's. The word "resorbable" means that the material is replaced during the period of implant deterioration by progressive reconstruction of the trabecular bone structures.

Resorbable implants have shown many advantages:

- make bone consolidation easier
- avoid the removal surgery
- avoid possible artifacts due to the presence of metal material during the MRI control
- avoid the metal ions release in bone tissues.

These implants are indicated in elective hand and foot surgery, traumas close to joints, pediatric surgery.

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Materials

RSB products are made of a compound (copolymer) of two lactide isomers: Poly (L-lactide-co- D,L lactide) more frequently called P(L/DL)LA 70:30 in literature, in conformity to ASTM F1925.

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The two isomers L-lactide and D-lactide are combined in order to form a single co-polymer that has a quite high resistance and a deterioration and reabsorption rate suitable for the most muscle-skeletal applications.


Ratio between the two components has been chosen for ensuring material's resistance to mechanical stresses during its lifetime and its complete reabsorption in the body.

The mechanical properties and the reabsorption times are represented in the following chart.

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**Phases of degradation
of amorphous
biodegradable
implants and tissue
reactions according to
H. Pistner**

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1. Healing phase: unchanged implant, development of a fibrous capsule with a high amount of fibroblasts.
2. Latent phase: unchanged implant, fibrous capsule gets thinner with less cells and more fibers or direct implant contact to bone structure.
3. Protracted resorptive phase: implant degradation begins with bulk erosion crack development, moderate cellular response with invasion of macrophages and foreign-body giant cells.
4. Progressive resorptive phase: progressive degradation of the implant associated with local immune response (macrophages, giant cells).
5. Recovery phase: no polymer remnants detectable, development of scar tissue and osseous replacement of the former implant site.

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Vite riassorbibile

CLINICAL INDICATIONS

Corrective osteotomies for hallux valgus according to Chevron-Austin surgical technique.

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


1.

Surgical technique

CHEVRON OSTEOTOMY OF THE FIRST METATARSUS

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- Establish the entry point for fixation at least 5 mm from the proximal osteotomy gap and introduce the guide wire Ø 1 mm HME 353-02 from dorsal to plantar up to the second cortex. Insert a temporary stabilization wire of Ø 1 mm, for fragment stabilization.

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- Check the accurate placement of the guide wire HME 353-02 and measure its depth with the HME 353-01 measuring device to determine the appropriate screw length.

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- Advance the guide wire through the distal cortex to minimize the risk of its accidental back sliding while drilling. Drill over the guide wire with the Ø 2.5 mm drill bit RSB 461-25 up to the previously determined screw length.

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
Surgical technique**Tecnica chirurgica**

4.

4. Assemble the tap for tip thread RSB 040035 to the straight cannulated handle STR 13A129 and ream distally over the guide wire up to the previously determined screw length, which can be read on the graduated scale on the tap stem.

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5.

5. Assemble the tap for head thread RSB 040045 to the straight cannulated handle STR 13A129 and ream over the guide wire up to the contact with the first cortex. Remove the guide wire HME 353-02.

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6.

6. Assemble the triangular screwdriver 1.5 mm RSB 030155 to the torque limiting handle STR 13D142 and screw the selected screw; if the screw protrudes dorsally, it is recommended to remove it temporarily, tap again with the tap for head thread and reinser the screw. Then remove the stabilization wire.

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Resorbable staple**Cambra riassorbibile****CLINICAL INDICATIONS**

Corrective osteotomies for hallux valgus according to Akin surgical technique and small fragments surgery.

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Resorbable pin**Pin riassorbibile****CLINICAL INDICATIONS**

Osteotomies and arthrodesis in foot and small fragments surgery.

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RSB 025-60

2,4 mm


Length 60 mm
Helicoidal profile for a better press-fit.

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

Surgical technique

STAPLE - EXAMPLE OF APPLICATION IN AKIN OSTEOTOMY OF THE FIRST PHALANX

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1. Place the drilling guide RSB 201 symmetrically across the osteotomy gap. Insert the drill bit Ø 2 mm 3600-020 in one of the 2 holes of the drilling guide and drill the first hole.

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2. Insert a Kirschner wire Ø 2 mm or the expulsor device RSB 205 through the drilling guide into the first drilled hole in order to keep the guide aligned. Insert the drill bit Ø 2 mm 3600-020 in the free hole and drill the second hole. Remove the Kirschner wire or the expulsor device and then the drilling guide.

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3. Decorticate up to the periosteum with a chisel to create the site of the staple bridge. Assemble the staple to the initial impactor RSB 202 and position the tips of the staple legs into the drilled holes. Advance the staple until the initial impactor is seated against the bone.

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4.

Surgical technique

4. Insert the expulsor device RSB 205 in the appropriate groove on the initial impactor and disassemble the staple from the latter.

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5.

5. Use the final impactor RSB 203 to fully place the staple against the cortical surface. If necessary disassemble the staple from the final impactor with the help of the expulsor device RSB 205 by inserting it in the appropriate groove.

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6.

6. For a greater osteotomy stability, a multiplanar insertion of a second staple is possible.

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Surgical technique

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PIN - EXAMPLE OF APPLICATION ON HAMMER TOE

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1.

- Resect the interphalangeal joint and prepare the pin site.

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2A.

- Insert the drill bit into the central portion of the proximal phalanx and drill first (fig. 2A), then into the central portion of the middle phalanx and drill beyond the distal phalanx (fig. 2B).

In case of normal or very hard bone it is recommended to use the drill bit with the same diameter of the pin to implant (drill bit 3600-018 for pin Ø 1.8 mm; drill bit 3600-025 for pin Ø 2.4 mm); in case of cartilage or cancellous bone, it is recommended to use the drill bit with the next smaller diameter than the one of the pin to implant (tip 3600-015 for pin Ø 1.8 mm; tip 3600-020 for pin Ø 2.4 mm).

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2B.

- Set the pin into the appropriate inserter RSB 00100 for the Ø 1.8 mm pin and RSB 00101 for the Ø 2.4 mm pin, taking care not to accidentally let it come out. Place the inserter on the middle phalanx and fully advance the pin by pushing the graduated bar until its end point.

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
5.

5. Place the inserter on the part of the pin that came out and push on the graduated bar to advance the pin from distal to proximal into the proximal phalanx hole, by inspecting visually the actual centering of the hole and its accurate insertion.

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6. Cut the excess part of the pin with an orthopaedic nipper and fully sink it in the bone structure by using the inserter. To avoid distal interphalangeal arthrodesis, it may be appropriate to break manually the pin by the aforementioned joint.

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Ghosh S, Viana JC, Reis RL, Mano JF.

Effect of processing conditions on morphology and mechanical properties of injection- molded poly(L-lactic acid) 29 May 2007 <https://doi.org/10.1002/pen.20799>.

Landes CA, Ballon A, Roth C.

In-patient versus in vitro degradation of P(L/DL)LA and PLGA J Biomed Mater Res B Appl Biomater. 2006 Feb;76(2):403-11.

RSB SCREW - VITE RSB

Barca F, Busa R.

Austin/chevron osteotomy fixed with bioabsorbable poly-L-lactic acid single screw J Foot Ankle Surg. 1997 Jan-Feb;36(1):15-20; discussion 79-80.

Wang Y, Song M, Xu Y, He X, Zhu Y.

Absorbable scaphoid screw development: a comparative study on biomechanics Ther Clin Risk Manag. 2016 Apr 22;12:643-50.

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Successful treatment of scaphoid fractures and nonunions using bioabsorbable screws: report of six cases

J Hand Surg Am. 2004 Jan;29(1):68-73.

RSB PIN - PIN RSB

Claes LE, Ignatius AA, Rehm KE, Scholz C.

New bioresorbable pin for the reduction of small bony fragments: design, mechanical properties and in vitro degradation Biomaterials. 1996 Aug;17(16):1621-6.

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Absorbable pins of self-reinforced poly-L-lactic acid for fixation of fractures and osteotomies J Bone Joint Surg Br. 1992 Nov;74(6):853-7.

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Chevron osteotomy of the first metatarsal stabilized with an absorbable pin: our 5-year experience Foot Ankle Int. 2013 Mar;34(3):380-5.

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Barca F, Busa R.

Resorbable Poly-L-Lactic Acid Mini-Staples for the Fixation of Akin Osteotomies. The Journal of Foot and Ankle Surgery 36(2):106-111 ,1997).

 Info

LSM-Med Srl
reserves the right
to make changes.

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**OVERMED Srl**