

LSMMed

Epibloc.1

**Product description
& surgical technique**

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
HME 0100 DISTAL RADIUS FRACTURES

A.O. classification:

- meta-epiphyseal fracture type A;
- type C (all) if the preliminary reduction is adequate.

Contraindications:

- fracture type B;
- fracture type C 3-3 with metaphyseal comminution.

HME 0101 PROXIMAL HUMERUS FRACTURES

Neer classification:

- 2-fragments fractures (anatomical or surgical neck with dislocation)
- 3-fragments
- 4-fragments if the preliminary reduction is adequate

Contraindications:

- 2-fragments (tuberosities fractures)
- 4-fragments anatomical neck with epiphyseal sinking or fragmentation

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Material

Wires: stainless steel AISI 316 L in conformity to ISO 5832-1

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


Pre-op view.

Controllo preoperatorio.

Post-op implant RX check and wrist joint mobility.

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Post-op view at 1 month and wrist joint mobility.

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
Technical features

- Epibloc.1 system uses the principle of the elastic endomedullary synthesis via percutaneous insertion of metallic wires in the medullary canal. The wires protruding out of skin are bended, strained and stabilized with an external plate.
- Available in 2 different versions: the one with Ø 2 mm wires, length 300 mm for distal radius fractures and the other with Ø 2.5 mm wires, length 500 mm for humerus fractures.
- The possibility to bend wires provides greater adaptability of the system to the specific clinical case in comparison to the use of the external fixators with fiches.
- The use of the polymeric plate reduces the dimensions of the external body compared to external fixators.
- Besides the 2 longitudinal grooves for the standard assembly, the plate has also 2 cross grooves to allow the possible insertion of a third and a fourth wire for the stabilization and the reduction of the most unstable fractures.

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Strumentario

Instrument
Strumentario


STR 13R127	Cannulated straight handle * Impugnatura retta cannulata *	
STR 031250	Cannulated screwdriver * Cacciavite cannulato *	Ø 2.5 mm

* 1/4" coupling / attacco 1/4"

∅ Hegaxon / Esagono

Surgical technique

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1. After fracture reduction, percutaneously insert the 2 wires inside the diaphysis and slide them pass along the walls of the two corticals opposite the point of entry (fig.1a, 1b).

2. Bend the protruding portions of the wires with appropriate forceps near the skin, till they are in opposite position, by performing the first bend at 90° and the second one, orthogonal to the first, at 100° (fig.2).

3. Tension the wires away from each other by distancing no more than 4 mm on each side; place them in the appropriate grooves of the plate and secure them permanently by tightening the plate locking screws with the appropriate screwdriver (Fig. 3). The correct tension makes the fracture sufficiently stable to allow immediate active mobilisation.

4. Perform a check of the correct reduction through the image intensifier. In case of wrong reduction, loosen the blocking screws of the plate, settle the wire force and reblock the system.

5. Cut off the portions of the wires exceeding the plate with appropriate cutters. Treat the skin in the entry points of the wires.

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Distal radius fractures: Treatment using the Epibloc™ system

Poggi D, Rocchi L, Merolli A, Fanfani F, Massarella M, and Catalano F.
Orthopaedics & Traumatology: Surgery & Research (2010) 96, 185-189.

Critical evaluation of 100 fractures of the wrist treated with the Epibloc System

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LSM-Med Srl
OVERMED Srl