

Product packaging

Packing materials common to all types of packaging		
Component of the packaging	Coding	Description
SHRINK FILM	PP 5	Plastic (Polypropylene (PPE))
BOX	PAP 21	Paper (Non corrugated cardboard)
LABEL	PET 1	Plastic (Polyethylene terephthalate (Polyester))
INSTRUCTIONS FOR USE	PAP 22	Paper

Sterile products packed in double blister		
Component of the packaging	Coding	Description
OUTER BLISTER	PETG 7	Plastic (PETG)
TYVEK OUTER SHEET	HDPE 2	Plastic (Tyvek ®)
PROTECTIVE SPONGES	LDPE 4	Plastic (Polyethylene)
INNER BLISTER	PETG 7	Plastic (PETG)
TYVEK INNER SHEET	HDPE 2	Plastic (Tyvek ®)

Sterile products packed in double vacuum pouch		
Component of the packaging	Coding	Description
OUTER POUCH	LDPE / PA 7	Plastic (OPA-PE)
INNER POUCH	LDPE / PA 7	Plastic (OPA-PE)
CONE COVER	OTHER 7	Plastic (Polycarbonate)

Sterile products packed in double Tyvek ® pouch		
Component of the packaging	Coding	Description
OUTER POUCH	HDPE / PET-PE	Plastic (Tyvek® / PET-PE film)
INNER POUCH	HDPE / PET-PE	Plastic (Tyvek® / PET-PE film)
PROTECTIVE POUCH	LDPE / PA 7	Plastic (OPA-PE)

Non-sterile products		
Component of the packaging	Coding	Description
POUCH	LDPE 4	Plastic (Polyethylene)
BLISTER	PETG 7	Plastic (PETG)
COMPOUND PACKAGING (Products with large dimensions)	C/LDPE 4	Pluriball – Adhesive tape
SEAL "NON STERILE"	PAP 22	Paper

Transport packaging

All the products		
Component of the packaging	Coding	Description
CARTON	PAP 20	Paper (Corrugated cardboard)
PLURIBALL	HDPE 2	Plastic (Polyethylene)
ADHESIVE TAPE	C/PVC 7	Plastic (PVC)
THERMAL CONTAINER	PS 6	Polistirolo